
Simpurity[™] Fibergel Ag Wound Dressing

Safe n Simple enhances quality of life with exceptional products at an affordable price


Introducing

Fibergel Ag with Gelling Action



Enhanced performance from a next generation silver fibergel dressing







Size	ltem#	Box/Case	HCPC
2"x2"	SNS56712	10 bxs/cs	A6196
4"x4"	SNS56744	10 bxs/cs	A6196
4"x4.75"	SNS56716	10 bxs/cs	A6197

KEY FEATURES:

- High silver content improved antimicrobial performance
- Fast initial silver elution quickly establishes an effective level of antimicrobial activity
- Sustained silver release on-going silver ion release allowing antimicrobial effects for up to 7 days
- High wet strength can be removed intact minimizing discomfort at dressing change

INDICATIONS:

For the management of moderate-to-highly exuding wounds such as:

- Pressure Ulcers
- Leg Ulcers
- Cavity Wounds
- Post operative or trauma wounds

For further information on this unique product please call sales at **1-844-767-6334**



Simpurity[™]

Safe n Simple enhances quality of life with exceptional products at an affordable price

HCPC

A6196

A6196

A6197

Fibergel Ag with Gelling Action

Enhanced performance from a next generation silver fibergel dressing

Size	ltem#	Box/Case
2"x2"	SNS56712	10 Pads/10 boxes
4"x4"	SNS56744	10 Pads/10 boxes
4"x4.75"	SNS56716	10 Pads/10 boxes

Product Details:



	Simpurity Fibergel Ag	Aquacel Ag Extra	Aquacel Ag
Material (Fibers)	Natural*	Chemically modified/processed	Chemically modified/processed
Silver content (excluding reinforcing)	2.4%	1.0%	1.0%
Silver Elution (at 12 hours)	49ppm	38ppm	30ppm
Remaining Silver Elution (at day 4)	40ppm	<2ppm	9ppm
Dressing class	Alginate	Hydrofiber	Hydrofiber
Haemostatic	Yes (minor bleeding)	No	No
Dry Strength (compared to Aquacel Ag Extra)	Good	Very good dry	Very inferior
Wet Strength (compared to Aquacel Ag Extra)	Intact on removal	Good weto	Very inferior
Gelling	Medium	High	High
Shrinkage	Nil	Up to 30% when wet	Up to 30% when wet
Felt weight g/100m ² (approx value)	155	175	120
Absorbency g/100cm² (approx value)	23	24	19

This data was prepared by an independent third party laboratory Simpurity[™] is a registered trademark of Safe n Simple, Ilc.

Aquacel® Ag and Aquacel® Ag Extra are registered trademarks of ConvaTec LTD.

SafenSimple WOUNDCARE 5827 Terex, Clarkston, MI 48346 | Toll Free: 844-767-6334 | www.sns-medical.com